


***ResMed***

## User Guide

English • Deutsch • Français • Italiano •  
Español • Português • Nederlands •  
Ελληνικά • Polski • Magyar • Русский •  
Slovenčina • Türkçe • Česky




**AirFit™** P10  
Nasal Pillows System

# AirFit™ P10


# AirFit™ P10

Fitting / Anlegen / Mise en place / Applicazione / Colocación / Colocação / Orzetten / Τοποθέτηση της μάσκας / Zakładanie maski / A maszk felhelyezése / Установка маски / Nastavenie masky / Maskenin Takılması / Nasazení masky

# AirFit™ P10

Adjustment / Einstellung / Ajustement / Regolazione / Ajuste / Ajuste / Afstellen / Προσαρμογή / Dopasowanie / Beállítás /  
Регулировка / Prispôsobenie / Ayarlama / Úprava nastavení masky



Removal / Abnehmen / Retrait / Rimozione / Para quitarla / Remoção / Afnemen / Αφαίρεση / Zdejmowanie / Levétel /  
Снятие маски / Snímanie / Çikarma / Sejmuti masky




# AirFit™ P10

Disassembly / Auseinanderbau / Démontage / Smontaggio / Desmontaje / Desmontagem / Demontage / Αποσυναρμολόγηση / Demontaž / Szétszerelés / Разборка / Rozobratie / Sökme / Demontáž

# AirFit™ P10

Reassembly / Zusammenbau / Remontage / Riassemblaggio / Montaje / Montagem / Opnieuw monteren / Επανα συναρμολόγηση / Składanie / Összeszerelés / Повторная сборка / Opätovné zloženie / Yeniden Montaj / Opřetovné sestavení

## English

### NASAL PILLOWS SYSTEM

This document provides the user instructions for:

- AirFit P10 (blue headgear) and AirFit P10 for Her (pink headgear) collectively referred to as AirFit P10.

### Intended use

The mask directs airflow noninvasively to a patient from a positive airway pressure device such as a continuous positive airway (CPAP) or bilevel system.

The mask is intended for:

- Patients (>30 kg) for whom positive airway pressure has been prescribed
- single patient re-use in the home environment and multipatient re-use in the hospital/institutional environment.

### Using the mask

- Refer to the Technical specifications section in this user guide for mask setting options.
- Follow the instructions provided by your physician or sleep therapist.
- To fit, adjust, remove, disassemble or reassemble your mask, follow the illustrations at the front of this guide.
- Do not overtighten the headgear straps. If you feel air leaks you may need to:
  - refit or reposition your mask
  - adjust your headgear
  - check that your mask is assembled correctly
  - contact your physician or sleep therapist to have your mask refitted.

- For a full list of compatible devices and mask setting options, see the “Mask/Device Compatibility List” on [www.resmed.com](http://www.resmed.com) on the **Products** page under **Service & Support**.
- The mask system does not contain natural rubber latex, PVC, DEHP or phthalates materials.

### Cleaning the mask at home

- The mask and headgear can only be handwashed.
- The mask vent may discolour over time. Follow the cleaning instructions in this guide to maintain performance and safety of the mask.

#### Daily/After each use

- To optimise the mask seal, remove facial oils from the pillows before use.
- Handwash the mask by rubbing gently in warm (30°C) water using mild soap.
- Rinse all components well with drinking quality water and allow them to air dry out of direct sunlight. Dry the vent by gently wiping with a cloth.
- If the vent is blocked, clean it with a soft bristle brush.
- Inspect each component and, if required, repeat washing until visually clean.
- When all components including the vent are dry, reassemble according to the reassembly instructions.

#### Weekly

Handwash the headgear in warm (30°C) water using mild soap. It may be washed without being disassembled. Rinse well and allow it to air dry out of direct sunlight.


## Reprocessing the mask between patients

The mask should be reprocessed when used between patients. Cleaning, disinfection and sterilisation instructions are available from the ResMed website, [www.resmed.com/masks/sterilization](http://www.resmed.com/masks/sterilization). If you do not have internet access, please contact your ResMed representative.

## Technical specifications

### Pressure-flow curve

The mask contains passive venting to protect against rebreathing.

AirFit P10 / AirFit P10 for Her	
Pressure	Flow (L/min)
4 cm H <sub>2</sub> O	20
8 cm H <sub>2</sub> O	29
12 cm H <sub>2</sub> O	37
16 cm H <sub>2</sub> O	43
20 cm H <sub>2</sub> O	49

### AirFit P10 / AirFit P10 for Her

**Deadspace:** Empty volume of mask to swivel. The deadspace of the mask varies according to pillows size.

123 mL (Large)

**Therapy pressure** (cm H<sub>2</sub>O)

4-20

**Resistance:** Drop in pressure measured (nominal) (cm H<sub>2</sub>O)

at 50 L/min 0.4

at 100 L/min 1.4

**Sound** – DECLARED DUAL-NUMBER NOISE EMISSION VALUES in accordance with ISO 4871. The A-weighted sound power level and the A-weighted sound pressure level of the mask, at a distance of 1 m, with uncertainty of 3 dBA, are shown.

Power level 21 dBA

Pressure level at 1 m 13 dBA

**Gross Dimensions** – Mask fully assembled – no headgear.

Height 590 mm

Width 149 mm

Depth 30 mm

**Environmental conditions** – all masks:

**Operating temperature:** +5°C to +40°C

**Operating humidity:** 15% to 95% non-condensing

**Storage and transport temperature:** –20°C to +60°C

**Storage and transport humidity:** up to 95% non-condensing

## Storage

Ensure that the mask is thoroughly clean and dry before storing it for any length of time. Store the mask in a dry place out of direct sunlight.

## Disposal

The mask does not contain any hazardous substances and may be disposed of with your normal household refuse.

## Symbols

The following symbols may appear on the mask or packaging.

 Caution, consult accompanying documents;

 Batch code;  Catalogue number;  European

Authorised Representative;  Temperature limitation;

 Humidity limitation;  Not made with natural rubber latex;

 Manufacturer;  Indicates a Warning or Caution and alerts

you to a possible injury or explains special measures for the safe and

effective use of the device;  Keep away from rain;  This way up;

 Fragile, handle with care;  Polypropylene recyclable.

## Warnings

- The mask must be used under qualified supervision for users who are unable to remove the mask by themselves. The mask may not be suitable for those predisposed to aspiration.
- The vent must be kept clear.
- The mask should only be used with CPAP or bilevel devices recommended by a physician or respiratory therapist.

- The mask should not be used unless the device is turned on. Once the mask is fitted, ensure the device is blowing air.
- Follow all precautions when using supplemental oxygen.
- Oxygen flow must be turned off when the CPAP or bilevel device is not operating so that unused oxygen does not accumulate within the device enclosure and create a risk of fire.
- Oxygen supports combustion. Oxygen must not be used while smoking or in the presence of an open flame. Only use oxygen in well ventilated rooms.
- At a fixed flow rate of supplemental oxygen flow, the inhaled oxygen concentration varies depending on the pressure settings, user breathing pattern, mask, point of application and leak rate.
- Your physician or sleep therapist is responsible for ensuring that the mask is compatible with the device. If used outside specification or if used with incompatible devices, the seal and comfort of the mask may not be effective and optimum therapy may not be achieved.
- Always follow cleaning instructions and use mild soap. Some cleaning products may damage the mask, its parts and their function, or leave harmful residual vapours that could be inhaled if not rinsed thoroughly.
- Do not combine components from different mask systems.
- Discontinue using the mask and consult your physician or sleep therapist if you have ANY adverse reaction.
- Using a mask may cause tooth, gum or jaw soreness or aggravate an existing dental condition. If symptoms occur, consult your physician or dentist.
- As with all masks, some rebreathing may occur at low CPAP pressures.
- Refer to your CPAP or bilevel device manual for details on its settings and operational information.

- Remove all packaging before using the mask.

### **⚠ Cautions**

- If any visible deterioration of a mask component is apparent (cracking, tears, etc), the mask component should be discarded and replaced.
- Avoid connecting flexible PVC products (eg, PVC tubing) directly to any part of the mask. Flexible PVC contains elements that can be detrimental to the materials of the mask, and may cause the components to crack or break.

### **Consumer warranty**

ResMed acknowledges all consumer rights granted under the EU Directive 1999/44/EG and the respective national laws within the EU for products sold within the European Union.

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© 2013 ResMed Ltd.
For patent information, see [www.resmed.com/ip](http://www.resmed.com/ip). AirFit is a trademark of ResMed Ltd.
**Manufacturer:** ResMed Ltd 1 Elizabeth Macarthur Drive Bella Vista NSW 2153 Australia.
**Distributed by:** ResMed Corp 9001 Spectrum Center Boulevard San Diego CA 92123 USA.
**EC REP** ResMed (UK) Ltd 96 Jubilee Ave Milton Park Abingdon Oxfordshire OX14 4RW UK.