



F5A Full Face Mask

Enlarge your vision by efficient face mask

F5A is one of the newest full face series. It is capable of auto adaptation to the nose bridge. Simple structure and light weight provide you with more comfortable and effective therapy.



Auto-adaptation

The groove on the cushion is able to adapt itself to the nose bridge, maintaining low air leakage and limited face pressure.

Liberty of vision

Based on the forehead frame-free design, users can embrace more vision with ease.

Anti-asphyxia valve

Under the premise of security, the position of anti-asphyxia valve is adjusted in order to achieve more functions on a smaller elbow.

More comfort

Due to the single-layer transitional design, the cushion touches the nose and face softly.

Product Code

F5A Full Face Mask

10030322 (L)
10030323 (M)
10030324 (S)



Replacement Headgear for F5A



Replacement Cushion for F5A

- Large size
- Medium size
- Small size

