

VisionAire™ 5

Stationary Oxygen Concentrator

CAIRE®
O₂ is what we do.



TIME-PROVEN DESIGN

PATIENT SATISFACTION

EASE OF SERVICE



The foundation of your patients' oxygen therapy

Compact and lightweight oxygen delivery

The VisionAire 5 has been completely redesigned to be state-of-the-art, meeting all current compliance standards with quality engineering and reliability.

Time-proven design

- **Industry-leading reliability**—greater than 99%* reliability rate out of the box
- **Comprehensive warranty**—3-year warranty on the device
- **Proven durability**—engineered with robust valves and sieve beds to minimize failure rates

Promotes patient satisfaction

- **Effective oxygen delivery**—up to 95.5% concentration from 1.0 to 5.0 LPM, meeting a wide variety of patient prescriptions
- **Quiet**—you will hardly know this lightweight and compact device is there
- **Energy efficient**—reduced power consumption to minimize the burden on the patient's energy bill

Ease of service

- **Maintenance free**—no required preventative maintenance reducing your total cost of ownership
- **Dual hour meter**—clear LCD provides time of use in hour and minutes
- **Enhanced diagnostics**—on board alarm alerts facilitate the ability to quickly troubleshoot



*Data on file.

VisionAire 5 Specifications

FLOW RATE

1.0 – 5.0 LPM

DIMENSIONS

52.8 x 35.8 x 29.2 cm

WEIGHT

13.6 kg

SOUND LEVEL

39 dB(A)

POWER

290 W (all flow rates)

MAXIMUM TUBING

15.2 m (50 ft)

WARRANTY

3 years



Accessories

NASAL CANNULA

Adult

7 ft. (2.1 m) CU002-1

25 ft. (7.6 m) CU002-4

Pediatric

7 ft. (2.1 m) CU002-7



HUMIDIFIER ADAPTOR

(9 in. / 23 cm.)

Item Number 20843882



HUMIDIFIER BOTTLE

Item Number HU003-1



OXYGEN TUBING

7 ft. (2.1 m) CU004-2

25 ft. (7.6 m) CU004-3

50 ft. (15.2 m) CU004-1

Serving your patients across the continuum of CAIRE.

Non-Delivery Solution

Reduce operating expenses and time in the field while improving patient care with the CAIRE non-delivery solution—the combination of the industry-leading FreeStyle® Comfort® portable oxygen concentrator with the maintenance-free VisionAire 5 stationary oxygen concentrator.



O₂ is what we do.

For more than 50 years, CAIRE has been designing and delivering the best technologies and services to support oxygen therapy.

CAIRE products are trusted globally by patients, providers, clinicians and the military. Let us put our experience to work for you.



Find out more about VisionAire 5 and accessory products at

CAIREinc.com

 CAIRE®

SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION. Please consult the applicable product instructions for use for product indications, contraindications, warnings, precautions, and detailed safety information.

© 2022 CAIRE Inc. All Rights Reserved. ML-CONC0082 G