




Simply CapFit™


Just like putting on a cap

Unlike any other mask, F&P Evora™ incorporates CapFit headgear that has been designed to be put on like a cap - a simple and intuitive movement that patients are familiar with.

The headgear structure ensures the seal is presented in the right place every time, making life easy for you and your patients.



Found Evora
simple to fit.¹




Found Evora simple to
take off and put on
in the dark.²




The next generation of Dynamic Support

The stability wings work in synergy with the floating seal to allow freedom of movement while keeping the mask comfortably in place.



Floating seal

Engages to wrap around the nose, providing a flexible and comfortable seal.



Stability wings

Designed to keep the floating seal in place throughout the night.



Supporting your business efficiency


Clinical research shows complexity in mask design³, poor mask fit, and discomfort⁴ can lead to more time spent with patients and ultimately non-adherence.⁴

When choosing a new mask consider how your choice can impact the amount of time you spend on:

- Fitting
- Educating
- Cleaning
- Ongoing support and motivation
- Resolving mask problems
- Subsequent mask setups



CapFit headgear
Minimal, one-size headgear for a simple and intuitive fit.



Soft-knit fabric
Comfortable headgear designed with AirEdges to avoid leaving marks on the face.



Compact nasal seal
Soft floating seal that sits comfortably under the nose to create an unobtrusive seal.



Radial air-flow
Laser-drilled exhaust holes designed to minimise noise and draft.



Fitting summary

Adjust the top strap



Adjust the back strap



Check the space

For further instructions, please refer to Fitting Your Mask section of the Evora Use and Care Guide

U Model

Code	Description
EVOISU	Evora Nasal Mask Small U Model
EVOIMU	Evora Nasal Mask Medium U Model
EVO1LU	Evora Nasal Mask Large U Model
EVO1WU	Evora Nasal Mask Wide U Model

Spare Parts

Code	Description
400EVO111	Evora Nasal Seal Small Spare
400EVO112	Evora Nasal Seal Medium Spare
400EVO113	Evora Nasal Seal Large Spare
400EVO114	Evora Nasal Seal Wide Spare
400EVO121	Evora Nasal Headgear Spare
400EVO122	Evora Nasal Backstrap Spare
400FPH341	F&P Swivel White
400EVO171	Evora Nasal Tube & Frame Spare
400EVO131	Evora Nasal Mask - Mask No H/G Small*
400EVO132	Evora Nasal Mask - Mask No H/G Medium*
400EVO133	Evora Nasal Mask - Mask No H/G Large*
400EVO134	Evora Nasal Mask - Mask No H/G Wide*

* Not available in all territories.

1. 20 out of 21 respiratory therapists rated the simplicity of the Evora Nasal mask as very simple or simple to fit. Internal validation study conducted with 21 respiratory therapists in United States of America 2019. 2. 43 out of 44 participants were able to put on and take off the Evora Nasal mask once fitted, in the dark. Internal validation study conducted with 44 participants in United States of America 2019. 3. Brown DL, Concannon M, Kaye AB, Zupancic M, Lisabeth LD. Comparison of two headgear systems for sleep apnea treatment of stroke patients. *Cerebrovasc Dis.* 2009;27(2):183-6. 4. Bachour A, Vitikainen P, Maasilta P. Rates of initial acceptance of PAP masks and outcomes of mask switching. *Sleep Breath.* 2016 May;20(2):733-8. 5. Weaver, T. E., & Grunstein, R. R. Adherence to continuous positive airway pressure therapy: The challenge to effective treatment. *Proc Am Thorac Soc.* 2008 Feb 15;5(2):173-8.

CapFit and Evora are trademarks of Fisher & Paykel Healthcare Limited. For patent information, please see www.fphcare.com/ip

REF 618459 REV B 2020-01 © 2020 Fisher & Paykel Healthcare Limited.